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- Ni, W. et al. The Serotonin Transporter is Present and Functional in Peripheral Arterial Smooth Muscle. *Journal of Cardiovascular Pharmacology*. 2004. 43: 770-781.
- B. Brenner, J.T. Hamey, B.A. Ahmed, B.C. Jeffus, R. Unal, J.L. Mehta, and F. Kilic. Plasma serotonin levels and the platelet serotonin transporter. *J Neurochem*. (2007) 102: 206-215.
- Carrasco G, Cruz MA, Gallardo V, Miguel P, Lagos M, and Gonzalez C (1998) Plasma and platelet concentration and platelet uptake of serotonin in normotensive and preeclamptic pregnant women. *Life Sci* 62: 1323-1332.
- Brenner B, Hamey JT, Ahmed BA, Jeffus BC, Unal R, Mehta JL, Kilic F (2007) Plasma serotonin levels and the platelet serotonin transporter. *J Neurochem* 102: 206-215.
- Watts SW and Fink G (1999) 5-HT<sub>2B</sub>-receptor antagonist LY-272015 is antihypertensive in DOCA-salt-hypertensive rats. *Am J Physiol Heart Circ Physiol* 276: 944-952.
- Watts SW (2005) 5-HT in systemic hypertension; foe, friend or fantasy? *Clin Sci* 108: 399-412.
- Amsterdam JD, Garcia-Espana F, Fawcett J, Quitkin FM, Reimherr FW, Rosenbaum JF, and Beasley C (1999) Blood pressure changes during short-term fluoxetine treatment. *J Clin Psychopharmacol* 19: 9-14.
- Diaz, J. Ni, W. Thompson, J, King A, Fink G, Watts SW. (2008) 5-hydroxytryptamine lowers blood pressure in normotensive and hypertensive rats. *JPET*. 325: 1031-1038

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2012.07 (rev01)

thoracic vena cava were harvested for basal 5-HT measurements.

#### iPRECIO® micro infusion pump:

During a second surgery using a similar anesthetic regimen, the iPRECIO® pump (Primetech, Tokyo Japan) was implanted through a small incision made at the base of the neck. Using blunt dissection a small subcutaneous pocket between the scapulae was created. (Figure 3) The pump was inserted and the skin sutured closed. Post recovery 5-HT or vehicle was loaded on each pump. The 5-HT solution was assessed by a two-way mixed design ANOVA and post-hoc testing at each time point with percreatinine complex and 1% ascorbic acid, antioxidant (0.02 g/pump). The solution was dissolved in 1N HCl, and a pH-balance (~7) is achieved with 4N NaOH.


Figure 3:

#### Results

The iPRECIO® micro infusion pump offered our laboratory an innovative tool to investigate the effects of delivering increasing doses of 5-HT to a single, conscious, freely moving, animal. We hypothesized that delivery of increasing doses of 5-HT to the S.D. rat would result in progressively greater reductions in BP when compared to the previous dose(s). This study demonstrated that increasing 5-HT did not progressively lower BP and that a smaller dose of 5-HT is required to elicit a hypotensive effect in the S.D. rats. (Figure 1& 2) These findings suggest that future experiments can be conducted at a lower dose, and is particularly relevant to studies focused on specific 5-HT agonists.

#### Discussion

This lab did discover that instead of a dose response to 5-HT, we were able to achieve an equivalent effect at a much lower dose than had been previously studied. This discovery would not have been possible in a fully ambulatory model with existing technology.

The iPRECIO® proved to be a technologic advancement from existing infusion methods and enabled this lab to make decisions moving forward that will save compound, labor and money while increasing physiologic relevance and overall animal welfare. This lab is pleased with the performance of the iPRECIO® pump and intends to use this product in further study.


Figure 2: . Effect of chronic 5-HT administration on mean arterial pressure (MAP) in 3 DOCA-salt hypertensive rats (red) and 5 sham normotensive rats (black).

1. Ni, W. et al. The Serotonin Transporter is Present and Functional in Peripheral Arterial Smooth Muscle. *Journal of Cardiovascular Pharmacology*. 2004. 43: 770-781.
2. Brenner, J.T. Harney, B.A. Ahmed, B.C. Jeffus, R. Unal, J.L. Mehta, and F. Kilic. Plasma serotonin levels and the platelet serotonin transporter. *J Neurochem*. (2007) 102: 206-215.
3. Carrasco G, Cruz MA, Gallardo V, Miguel P, Lagos M, and Gonzalez C (1998) Plasma and platelet concentration and platelet uptake of serotonin in normotensive and preeclamptic pregnancies. *Life Sci* 62: 1323-1332.
4. Brenner B, Harney JT, Ahmed BA, Jeffus BC, Unal R, Mehta JL, Kilic F (2007) Plasma serotonin levels and the platelet serotonin transporter. *J Neurochem* 102: 206-215.
5. Watts SW and Fink G (1999) 5-HT<sub>2B</sub>-receptor antagonist LY-272015 is antihypertensive in DOCA-salt-hypertensive rats. *Am J Physiol Heart Circ Physiol* 276: 944-952.
6. Watts SW (2005) 5-HT in systemic hypertension; foe, friend or fantasy? *Clin Sci* 108: 399-412.
7. Amsterdam JD, Garcia-Espana F, Fawcett J, Quitkin FM, Reimherr FW, Rosenbaum JF, and Beasley C (1999) Blood pressure changes during short-term fluoxetine treatment. *J Clin Psychopharm* 19: 9-14.
8. Diaz, J. Ni, W. Thompson, J. King A, Fink G, Watts SW. (2008) 5-hydroxytryptamine lowers blood pressure in normotensive and hypertensive rats. *JPET*. 325: 1031-1038